

SAMPLE PREPARATION AND SHIPPING INSTRUCTIONS

SEROLOGY

1. If sending serum samples that are not heat inactivated, ship in an insulated container on a cold pack. If sending serum samples that are heat inactivated, ship at ambient temperature in an insulated container.
2. Vials should be tightly sealed and labeled with ID number using permanent marker. Outside threaded plastic tubes preferred to minimize leakage.
3. Submission form and a purchase order must accompany each shipment.
4. Receiving Department is closed on Sundays and holidays.
5. Advance notification of all shipments allows for the most efficient processing of samples.

Ship Samples to:

Attn: VRL
401 Professional Drive
Suite 210
Gaithersburg, MD 20879

Questions: Call Client Services at 301-610-2521 or 800-804-3586, fax 240-686-6776 or visit www.vrl.net

Serology Requirements:

Sample Type	Amount Needed	Sample Requirements
Serum for Mouse Core Serology Panel	*Undiluted: 0.05 ml needed per profile Diluted: 0.25 ml needed per profile	Dilute 1 part serum/plasma in 4 parts sterile Phosphate Buffer Saline (PBS) Heat Inactivate @ 56°C
Rodent Serum	*Undiluted: 0.025 ml needed per assay 0.20 ml needed per profile Diluted: 0.05 ml needed per assay 1.0 ml needed per profile	Dilute 1 part serum/plasma in 4 parts sterile Phosphate Buffer Saline (PBS) Heat Inactivate @ 56°C
Rabbit Serum	*Undiluted: 0.025 ml needed per assay 0.20 ml needed per profile	Do not Heat Inactivate @ 56°C

